
NUCLEODUR® high purity silica for HPLC

NUCLEODUR® is a fully synthetic type B silica (silica of 3rd generation) offering highly advanced physical properties like **totally spherical** particle shape, outstanding **surface microstructure**, high **pressure stability** and **low metal content**.

NUCLEODUR® as a state-of-the-art silica is the ideal base material for modern HPLC phases. It is the result of MACHEREY-NAGEL's pioneering research in chromatography for more than 40 years and succeeds MN's famous NUCLEOSIL® silica.

In RP liquid chromatography the efficiency of the packing is strongly affected by the quality of the base silica itself. Shortcomings in the surface geometry of the particles or metal contaminants are the main reasons for inadequate coverage with the covalently bonded alkylsilanes in the subsequent derivatization steps. It is well known, that poor surface coverage and, in consequence, high activity of residual free silanols often results in peak tailing or adsorption, particularly with basic compounds.

Particle shape and surface symmetry

NUCLEODUR® silicas are synthesized in a unique and carefully controlled manufacturing process which provides silica particles, which are totally spherical. The picture shows the outstanding smoothness of the NUCLEODUR® surface.

Purity

As already mentioned above, a highly pure silica is required for achieving symmetric peak shapes and maximum resolution. Inclusions of e.g. iron or alkaline earth metal ions on the silica surface are largely responsible for the unwanted interactions with ionizable analytes, e.g. amines or phenolic compounds.

NUCLEODUR® is virtually free of metal impurities and low acidic surface silanols. Elemental analysis data of NUCLEODUR® 5 µm measured by AAS are listed below.

Elementary analysis (metal ions) of NUCLEODUR® 100-5

Aluminium	< 5	ppm
Iron	< 5	ppm
Sodium	< 5	ppm
Calcium	< 10	ppm
Titanium	< 1	ppm
Zirconium	< 1	ppm
Arsenic	< 0.5	ppm
Mercury	< 0.05	ppm

Pressure stability

The totally spherical and 100% synthetic silica gel exhibits an outstanding mechanical stability, even at high pressures up to 800 bar and elevated eluent flow rates.

In addition, after several cycles of repeated packing, no significant drop in pressure can be observed. The latter is of prime importance for preparative and process-scale applications.

Physical properties of NUCLEODUR®

Surface (BET)	340 m ² /g
Pore size	110 Å
Pore volume	0.9 ml/g

NUCLEODUR® modifications

Several different surface modifications based on NUCLEODUR® silica have been developed over the last years providing a full range of specified HPLC phases and an ideal tool for every separation:

- ✦ NUCLEODUR® C₁₈ Gravity and C₈ Gravity
- ✦ NUCLEODUR® C₁₈ Isis
- ✦ NUCLEODUR® C₁₈ Pyramid
- ✦ NUCLEODUR® Sphinx RP
- ✦ NUCLEODUR® CN and CN-RP
- ✦ NUCLEODUR® NH₂ and NH₂-RP
- ✦ NUCLEODUR® C₁₈ ec and C₈ ec

For important properties of NUCLEODUR® phases please see our summary.

Overview of NUCLEODUR® HPLC phases

Columns for HPLC

Phase	Specification	Characteristics*			Stability	Structure
		A	B	C		
C₁₈ Gravity	octadecyl phase, high density coating multi-encapping 18 % C · USP L1				pH stability 1 - 11, suited for LC/MS	NUCLEODUR® (Si-O ₂) _n
C₈ Gravity	octyl phase, high density coating multi-encapping 11 % C · USP L7				pH stability 1 - 11, suited for LC/MS	NUCLEODUR® (Si-O ₂) _n
C₁₈ Isis	octadecyl phase with specially crosslinked surface modification endcapping 20 % C · USP L1				pH stability 1 - 10, suited for LC/MS	NUCLEODUR® (Si-O ₂) _n
C₁₈ Pyramid	C ₁₈ modification with polar endcapping 14 % C · USP L1				stable in 100 % aqueous eluents without phase collapse, pH stability 1 - 9, suited for LC/MS	NUCLEODUR® (Si-O ₂) _n
Sphinx RP	bifunctional RP phase, balanced ratio of propyl-phenyl and C ₁₈ ligands; endcapping 15 % C; USP L1 and L11				pH stability 1 - 10, suited for LC/MS	NUCLEODUR® (Si-O ₂) _n
C₁₈ ec	octadecyl phase, medium density coating endcapping 17.5 % C · USP L1				pH stability 1 - 9	NUCLEODUR® (Si-O ₂) _n
C₈ ec	octyl phase, medium density coating endcapping 10.5 % C · USP L7				pH stability 1 - 9	NUCLEODUR® (Si-O ₂) _n
CN / CN-RP	cyano (nitrile) phase for NP and RP separations 7 % C · USP L10			-	pH stability 1 - 8, suited for mobile phases with high contents of water	NUCLEODUR® (Si-O ₂) _n
NH₂ / NH₂-RP	amino phase for NP and RP separations 2.5 % C · USP L8			-	pH stability 2 - 8, suited for mobile phases with high contents of water	NUCLEODUR® (Si-O ₂) _n
SiOH	unmodified USP L3	-	n.a.	-	pH stability 2 - 8	(Si-O ₂) _n

* A = hydrophobic selectivity, B = polar / ionic selectivity, C = steric selectivity

An optimised phase for every separations

Application	Similar phases**	Separation principle · Retention mechanism
in general compounds with ionizable functional groups such as basic pharmaceuticals and pesticides	NUCLEOSIL® C₁₈ HD Waters Xterra® RP ₁₈ / MS C ₁₈ ; Phenomenex Luna® C18 (2), Synergi™ and Max RP; Zorbax® Extend C18; Inertsil® ODS III; Purospher® RP-18, Star RP-18	only hydrophobic interactions (van der Waals interactions)
like C ₁₈ Gravity, however generally shorter retention times for nonpolar compounds	NUCLEOSIL® C₈ HD Waters Xterra® RP ₈ / MS C ₈ ; Phenomenex Luna® C8; Zorbax® Eclipse; XDB-C8	
high steric selectivity, thus suited for separation of positional and structural isomers, planar / non-planar molecules	NUCLEOSIL® C₁₈ AB Inertsil® ODS-P; YMC® Pro C18RS	steric interactions and hydrophobic interactions
basic pharmaceutical ingredients, very polar compounds, organic acids	Phenomenex Aqua®; YMC® AQ; Waters Atlantis® dC18	hydrophobic interactions and polar interactions (H bonds)
compounds with aromatic and multiple bond systems	no similar phases	π-π interactions and hydrophobic interactions
robust C ₁₈ phase for routine analyses	NUCLEOSIL® C₁₈ Spherisorb® ODS II; Hypersil® ODS; Waters Symmetry® C18; Inertsil® ODS II; Kromasil® C18; LiChrospher® RP 18	only hydrophobic interactions (van der Waals interactions)
robust C ₈ phase for routine analyses	NUCLEOSIL® C₈ ec / C₈ Spherisorb® C8; Hypersil® MOS; Waters Symmetry® C8; Kromasil® C8; LiChrospher® RP 8	some residual silanol interactions
polar organic compounds (basic drugs, molecules containing π electron systems)	NUCLEOSIL® CN / CN-RP	π-π interactions, polar interactions (H bonds), hydrophobic interactions
sugars, sugar alcohols and other hydroxy compounds, DNA bases, polar compounds in general	NUCLEOSIL® NH₂ / NH₂-RP	polar / ionic interactions, hydrophobic interactions
polar organic compounds in general	unmodified NUCLEOSIL®	polar / ionic interactions

** phases which provide a similar selectivity based on chemical and physical properties


Columns for HPLC

Particle size and separation efficiency

1.8 μm particles for increased separation efficiency

- decrease of analysis time (ultra fast HPLC)
- shorter columns with high separation efficiency
- significant improvement of resolution
- increased detection sensitivity
- suitable for LC/MS due to low bleeding characteristics
- all NUCLEODUR® premium phases are available in 1.8 μm : C₁₈ Gravity, C₈ Gravity, C₁₈ Isis, C₁₈ Pyramid, Sphinx RP
- NUCLEODUR® 1.8 μm particles are fractionated to limit the increase in back pressure

Now available: 1.8 μm particle size!

Features of 1.8 μm NUCLEODUR® silica particles

- increase of separation efficiency by higher number of theoretical plates (N)
- significant improvement in resolution
- low column back pressure

Comparison of back pressure:

Eluent: 100 % methanol
Flow rate: 1.5 ml/min
Temperature: 22 °C
Column dimension: 50 x 4.6 mm


	NUCLEODUR® C ₁₈ Gravity	Competitor A
3 μm	70 bar	-
1.8 μm	130 bar	170 bar

- shorter run times

Resolution as a function of particle size

Column: 50 x 4 mm NUCLEODUR® C₁₈ Gravity
A) 3 μm , B) 1.8 μm
Eluent: acetonitrile – water (80:20, v/v)
Flow rate: 2 ml/min
Pressure: A) 80 bar, B) 160 bar
Detection: UV, 254 nm


Peaks:
1. Naphthalene
2. Ethylbenzene

Columns for HPLC

Reduction of analysis time

Column: 50 x 4 mm (for 5 μm 125 x 4 mm) NUCLEODUR® C₁₈ Isis
Eluent: 100 % methanol
Flow rates and pressure see figure
Detection: UV, 254 nm

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Surface modification

By use of specific C₁₈ silanes and appropriate polymeric bonding technologies a dense shield of alkyl chains protects the subjacent silica matrix. Elemental analysis of NUCLEODUR® C₁₈ Isis shows a carbon load of 20%.

The target crosslinking of the C₁₈ chains on the surface enables the separation of compounds with similar molecular structure but different stereochemical properties. The technical term for this feature is steric selectivity.

The separation of *o*-terphenyl and triphenylene is a concrete example to evaluate the selectivity potential of a reversed phase column in terms of the different shape of two molecules. The phenyl rings of *o*-terphenyl are twisted out of plane while triphenylene has a planar geometry.

The separation factor (α value) is a measure for the steric selectivity. As is shown in the following chromatograms the α value is considerable larger on NUCLEODUR® C₁₈ Isis compared to a conventional C₁₈ column.

Steric selectivity of NUCLEODUR® C₁₈ Isis

Columns: 125 x 4 mm; **NUCLEODUR® C₁₈ Isis, monomerically coated C₁₈ phase, C₁₈ phase with polar endcapping**

Eluent: methanol – water (90:10, v/v)

Flow rate: 1 ml/min, temperature 35 °C

Detection: UV, 254 nm

Injection volume: 5 µl


Peaks:

1. *o*-Terphenyl
2. *m*-Terphenyl
3. *p*-Terphenyl
4. Triphenylene

C₁₈, polar endcapping

monomeric C₁₈

C₁₈ Isis

Steric selectivity of NUCLEODUR® C₁₈ Isis

Columns 125 x 4 mm; eluent methanol – water (80:20, v/v)


Flow rate: 1 ml/min, temperature 40 °C

Detection: UV, 254 nm, injection volume 1 µl

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MACHERY-NAGEL

RP HPLC with highly aqueous eluents

Conventional reversed phase columns often display stability problems in eluent systems with high percentage of water (> 95%) as evidenced by a sudden decrease of retention time and overall poor reproducibility. This phenomenon is described as phase collapse caused by the mobile phase expelled from the pores due to the fact, that hydrophobic RP phases are incompletely wetted with the mobile phase.

Different approaches can be used to increase column stability with highly aqueous mobile phase systems. The most promising concepts are incorporating a polar group in the hydrophobic alkyl chain, or using hydrophilic endcapping procedures to improve the wettability of the reversed phase modification.


Stability features

NUCLEODUR® C₁₈ Pyramid is a silica phase with hydrophilic endcapping, designed especially for use in eluent systems of up to 100% water. The stability test shows the retention behaviour of tartaric, acetic and maleic acid under purely aqueous conditions on NUCLEODUR® C₁₈ Pyramid in comparison with a conventionally bonded RP phase.

It can be shown that the retention times for NUCLEODUR® C₁₈ Pyramid remain nearly unchanged between initial injection and restart after the flow has been stopped for 12 hours, whilst the performance of the conventional RP column collapsed totally after 5 min.

Retention characteristics


The polar surface derivatization exhibits retention characteristics, which differentiate the "Pyramid" from conventional C₁₈ stationary phases. The chromatogram at right shows the improved retention behaviour of very polar compounds such as short chain organic acids, which are insufficiently retained on RP columns with predominantly hydrophobic surface properties.

Separation of very polar compounds

Column: 125 x 4 mm NUCLEODUR® C₁₈ Pyramid, 5 µm
Eluent: 0.2% H₃PO₄
Flow rate: 1,0 ml/min
Temperature: 22 °C
Detection: UV, 202 nm
Injection volume: 2 µl

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Alternative RP selectivity

NUCLEODUR® Sphinx RP is characterized by exceptional selectivity features generated by a well-balanced ratio of covalently bonded octadecyl and phenyl groups. The combination of classical hydrophobic with π - π interactions (aromatic ring system) expands the scope of selectivity in comparison with conventional reversed phase packings. NUCLEODUR® Sphinx RP is particularly suited for the separation of molecules containing aromatic and multiple bonds. For the separation of polar compounds NUCLEODUR® Sphinx RP can be especially recommended and can also outperform many customary C₁₈ phases.

In addition, exhaustive endcapping steps minimize unwanted surface silanol activity and guarantee excellent peak shapes even for strongly basic analytes.

Different from standard phenyl phases, NUCLEODUR® Sphinx RP is far more stable towards hydrolysis and is also suggested for LC/MS applications.

Due to the additional intermolecular interactions NUCLEODUR® Sphinx RP is an interesting replenishment to the high density bonded phases NUCLEODUR® C₈/C₁₈ Gravity and the polar endcapped NUCLEODUR® C₁₈ Pyramid.

Separation of flavonoids on 3 different NUCLEODUR® phases

Columns: 150 x 4.6 mm

A) NUCLEODUR® C₈ Gravity, 5 µm

B) NUCLEODUR® C₁₈ Gravity, 5 µm

C) NUCLEODUR® Sphinx RP, 5 µm


Eluent: water – methanol (40:60, v/v), 1 ml/min, 30 °C

Detection: UV, 270 nm; injection volume: 3 µl

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NUCLEODUR® C₁₈ ec for daily routine analysis and up-scaling in preparative HPLC

The efficiency of a separation is controlled by particle size and selectivity of the stationary phase. The exceptional surface coverage of monomeric bonded alkylsilanes, combined with an exhaustive endcapping, results in a surface with lowest silanol activity. This allows the tailing-free elution of polar compounds such as basic drugs. NUCLEODUR® C₁₈ ec is also ideal for scale-up purposes.

Chemical stability

The utmost purity of the base silica and the exceptional silane bonding chemistry minimizes the risk of dissolution, or hydrolysis at pH extremes.

High loadability

Loadability, probably the most important feature for preparative LC, is determined by pore size, pore volume and surface area of the packing.

NUCLEODUR® octyl phases

Based on the same totally spherical and highly pure silica the C₈ phases exhibit the same excellent chemical and mechanical stability features as the C₁₈ counterparts. Due to the shorter chain and less hydrophobic properties of the stationary phase the retention of nonpolar compounds is decreased, and in consequence a reduction in time of analysis can be achieved. Moreover a stronger polar selectivity, particularly with the separation of ionizable analytes is frequently observed (as distinct from the C₁₈ phases).


Some general principles are:

- High density C₈ and C₁₈ phases allow tailing-free elution even for very polar compounds
- Octyl phases (C₈) show superior polar selectivity
- Octadecyl phases (C₁₈) show superior hydrophobic selectivity
- Hydrophobic compounds show shorter retention times on C₈ phases

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
Available standard dimensions of VarioPrep columns with axially adjustable end fitting

ID [mm]	Length [mm]									End fitting design
	10*	15*	50	75	100	125	150	250	500	
8	X		X		X	X	X	X		
10			X		X	X	X	X		
16	X		X		X	X	X	X		
21			X	X	X	X	X	X		
32		X			X		X	X		
40			X		X	X	X	X	X	
50		X			X		X	X		
80							X	X		

* 10 x 8, 10 x 16, 15 x 32 and 15 x 50 mm ID columns are used as guard columns and require adequate holders.

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MACHEREY-NAGEL

Available standard dimensions of EC columns · please ask for availability of certain phases


ID [mm]	Length [mm]											End fitting design	
	8*	20	30	50	75	100	125	150	200	250	300		
2	-	X	X	X	X	X	X	X	X	X	X	X	
3	X	X	X	X	X	X	X	X	X	X	X		
4	X	X	X	X	X	X	X	X	X	X	X		
4,6	-	X	X	X	X	X	X	X	X	X	X		

* Please note that 3 mm ID guard column cartridges are applicable for 2 mm ID and 3 mm ID EC columns, while 4 mm ID guard column cartridges are also used for 4.6 mm ID EC columns.

Installation of the EC guard column adaptor (Cat.No. 7.081 898)

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