



Reactor-Ready Lab Reactor







Replace multiple reaction set-ups with a single universal system that accepts vessels from 100 ml to 5 litres









easy drain down
reduces stress on vessel sidearms

Change to Reactor-Ready


It improves your chemistry by saving space, time and money

Reactor-Ready Duo



- slide and glide
- set and lock





Patented stirrer coupling

- drop in
- no tools
- no fuss

Standard vessels

- 250 ml to 5 litres
- · vacuum jacketed option



Optional AVA Control Software allows you to log and control stirrers, circulators, balances, pumps,

temperature sensors



- mimic plant scale geometry
 - 100 ml to 5 litres
 - · vacuum jacketed option



Accepts all leading brands of overhead stirrer





Compatible with all leading brands of circulator or thermoregulator



Reactor-Ready

Swap reaction vessels in minutes, NOT hours

Reactor-Ready is designed as a universal reactor work station that can be used for different vessel sizes and different experiments. Reactor-Ready can easily and

quickly be configured to suit the chemistry and scale you need for each project.

The beauty of Reactor-Ready is that this one reactor work station can replace many, saving you money and fumehood space.

Features include

- Rapid vessel exchange with quick-release vessel clamp and hose couplings.
- Range of single and vacuum jacketed vessels from 100 ml to 5 litres.
- Range of process vessels to mimic larger scale plant or manufacturing reactors.
- · Accepts all leading brands of overhead stirrer and allows easy, tool-free adjustment.
- · Triple support stand features heavy duty stainless steel support rods for strength and stability.
- · Self aligning stirrer coupling engages without the need for tools.
- Temperature range: -70 °C to +230 °C (jacket temperature)
- Innovative hose manifolds allow easy thermofluid drain down.
- · Wide range of accessory glassware including condensers, dropping funnels etc.
- Optional software allows you to log and control stirrers, circulators, balances, pumps, temperature sensors and other devices.

Patented vessel clamp allows rapid vessel exchange without tools

- Unique clamp design allows vessel to be supported and clamped by a single operator in less than a minute.
- · Use the same clamp to hold all vessel sizes; plus the vessels remain supported even when the clamp is unlocked.













Stainless steel vessel support collar

Vessel support collar aids vessel loading

- Reactor-Ready uses an innovative support collar fitted to the neck of each vessel, which allows the vessels to be swapped in and out of the Reactor-Ready core like a cartridge.
- The stainless steel collar locks securely into the support plate on the framework and allows the vessel to remain supported, even when not clamped.
- Vessel exchange is tool free and can easily be undertaken by one person.

Unique self aligning stirrer coupling saves time when changing vessels

- Patented drop-in stirrer coupling engages easily without the need for tools.
- Simply lower the stirrer coupling to engage with the stirrer guide below.
- The stainless steel stirrer coupling will fit in any standard stirrer chuck (for 10mm diameter)
- Powerful stirring up to 500 rpm, or 800 rpm for short periods (subject to motor selected).









Choice of PTFE stirring shafts

- · Reactor-Ready is offered with a choice of PTFE stirring shafts: anchor, retreat curve and turbine/propeller styles.
- PTFE anchor stirring shaft is included with all standard reaction vessel kits
- PTFE turbine stirring shaft is included with all process reaction vessel kits.

I-beam support allows the stirring motor to be positioned with ease

- I-beam allows the stirrer motor to be raised or lowered without needing to remount the stirrer.
- Adjustable safety stop collars located below the I-beam allow the user to lower the stirrer motor to a set position for ease of operation and safety.
- · Reactor-Ready accepts all leading brands of motor or air-powered overhead stirrer with support rods up to 16 mm in diameter.





PTFE Stirrer Guide

- The Reactor-Ready's PTFE stirrer guide screws directly on to the glass lid and features a unique extended support which minimises wobble even at high speeds.
- The stirrer guide is chemically inert, gas-tight and suitable for work under vacuum.

PTFE stirrer guide provides whip-free agitation

PTFE anchor

stirrer



PTFF retreat

curve stirrer

Quick release PEEK hose couplings

- Operates over the temperature range -70 °C to +230 °C (jacket temperature).
- · Quick-release, screwthread design allows rapid exchange of vessels without the use of tools.
- · Simply tighten by hand for a leak-tight seal.
- PEEK with Viton O-rings ensures all wetted parts are chemically resistant.
- · Lightweight, metal-free design prevents chipping and minimises strain on the glass fitting.
- All hoses have a 12 mm ID, maximising thermofluid flow and heating/cooling performance.
- · Low thermal conductivity of PEEK reduces frosting when cooling.





PTFE turbine/

propeller stirrer

Hose manifolds allow easy thermofluid drain down

- The inlet and outlet hoses from the Reactor-Ready vessel connect to manifolds attached to the support rod. These manifolds serve two important functions:
 - 1. Support the hoses running from the vessel, thereby eliminating stress on the glass sidearms.
 - 2. Allows thermofluid drain-down, which provides a quick and clean method of removing the thermofluid from the vessel jacket when changing vessels.
- · By opening the venting valve on the (upper) outlet manifold, thermofluid drains from the vessel back to the thermoregulator under gravity.
- Opening the drain tap on the (lower) inlet manifold allows drainage to an external reservoir.





Case study: Accelerating process optimisation in a busy chemical development lab

Researchers at a leading US pharmaceutical company invested in Reactor-Ready for use in their chemical development lab. The process R&D group, specialising in small molecule research, are using Reactor-Ready for optimising new or existing chemistry, identifying and eliminating impurities, and screening process variables for transfer to kilo-labs, pilot plant and manufacturing sites.

Commenting on their existing methods the Lead Scientist said, "We use round bottom flasks, but temperature control via mantles or oil baths is poor. Alternatively, when we use traditional jacketed reactor vessels, we need a different set up for every scale change. Even changing a reactor of the same size takes ages."

"So I liked Reactor-Ready the minute I saw how easy it was to change reactors. Reactor-Ready does it all. I get one stand which fits all sizes up to 5 litres. I get overhead stirring on all scales. I get vacuum jacketed vessels, which provide excellent temperature control with minimum heat loss, which is a real advantage at low temperatures, and I get the ability to rapidly change out my reactors."

"The first time I changed vessels, it took about 5 minutes but now I can change in less than 2 minutes. That includes removing and attaching the heating fluid hoses. Occasionally I lose a drop or two of heat transfer fluid but most times I lose none. The thing I like about the system is the control and speed. I get better control of my reactions and I can vary my work or scale immediately".







Select the overhead stirrer you need.



4. If you need a thermoregulator, hoses, hose adapters, thermofluid or accessory glassware, then select from the accessory list.

5. If you need automation add the AVA Control Software and data hub.



Reactor-Ready in detail

It is the attention to detail and wealth of innovative features that makes Reactor-Ready so much better



Precision ground lid and vessel flanges combine with PTFE collar and FEP-encapsulated O-ring to offer superior sealing performance.

Vessel clamp support allows easy vertical repositioning of the vessel whilst maintaining the stirrer coupling alignment.



Venting stopcock on the outlet manifold allows quick and easy draining of the thermofluid during vessel exchange.

Vessel-to-manifold hoses are supported by the manifolds and frame, eliminating stress on the vessel.

> Quick-release PEEK coupling allows rapid exchange of vessels without the use of tools. Operates over the temperature range -70 °C to +230 °C (jacket temperature).






Thermofluid may be drained back to the thermoregulator or decanted via a fluid drain tap on the inlet manifold.


Heavy-duty base plate with levelling feet ensures that the Reactor-Ready support stand is secure and stable.






Reactor-Ready core system includes

- Stand with base, support rods, bosses and all clamps
- · Glass reaction vessel lid
- · Vessel FEP O-ring seal and PTFE collar
- HDPE drip tray
- · Vessel-to-manifold insulated hoses with quick-release connections
- Lower inlet manifold + drain





· Two stainless steel vessel support collars

· Locking rod for support rod bosses



Reaction vessel kits make ordering vessels simple and more cost effective

These cost effective vessel kits include the vessel and accessories you need to get started. The Pt100 probe and PTFE stirrer included in each kit are matched to the vessel size. Each vessel kit includes:

- · Reaction vessel with PTFE stopcock
- · Sidearm couplings Pt100 PTFE temperature probe with LEMO

connector PTFE temperature





Reactor-Ready Duo

Supports two glass reaction vessels of different or same volume from 100 ml to 5 litres

Reactor-Ready Duo shares the same unique features as Reactor-Ready, except that it supports two separate jacketed glass reaction vessels. The system can be configured to operate with a single thermoregulator controlling the temperature of both vessels simultaneously, or with two thermoregulators controlling the temperature of each vessel independently.



Unique vessel clampsupports two vessels

- change vessels independently
 - change vessels in minutes

Vessels

- 100 ml to 5 litres
- combine different volumes
- vacuum jacketed option
- process and standard geometry

Thermofluid manifolds

- choice of two manifold kits
- single circulator manifolds
- double circulator manifolds



Strong framework

- heavy duty base
- solid s/steel support rods
 - large drip tray

Why Reactor-Ready Duo?

All the benefits of Reactor-Ready with two vessels in parallel or series

Features include

- · Rapid exchange of both vessels independently, with quick-release vessel clamp and hose couplings.
- · Choice of manifold kits allow two vessels to run from a single thermoregulator or two separate thermoregulators.
- System accepts two overhead stirrers which can be operated and moved independently.

Applications

- Parallel synthesis or reaction optimisation: use the same or different size vessels and vary stirring speed, stirrer shape and temperature between vessels as appropriate.
- Two stage reaction: transfer reactant from one vessel to the other using vacuum or a pump.
- Single reaction vessel: using the second vessel as either a receiving or feed vessel (where reagents can be pre-heated or pre-cooled prior to addition).
- Use optional AVA Software to control fluid transfer between vessels.







Different vessels with










How to order Reactor-Ready Duo...







4. Select the overhead stirrers.



5. If you need a thermoregulator, hoses, hose adapters, thermofluid or accessory glassware, then select from the accessory list.

6. If you need automation, add the AVA Software and data hub.



Reactor-Ready Duo core system

Easy-to-order Reactor-Ready Duo core system includes all the framework components required for operation

Features include

 Robust heavy-duty framework features solid stainless steel support rods.

 Two separate quick release clamps support each of the vessels independently.

· Large drip tray will retain vessel contents.

· Manifold and hose kits must be ordered separately.

Core system (RR121300) includes

• Stand with base, support rods, bosses and all clamps

• 2 x glass reaction vessel lids

2 x vessel FEP O-ring seals and PTFE centering collars

• 2 x PTFE stirrer guides and couplings

- HDPE drip tray
- 2 x stainless steel stirrer drop-in couplings
- 2 x stainless steel vessel support collars
- locking rod for support rod bosses



Reactor-Ready Duo manifolds and hoses

A choice of manifold kits for use with one or two thermoregulators

Double thermoregulator manifold and hose kit (RR121310)
This double thermoregulator kit is designed to work with two
separate thermoregulators controlling the temperature of each

vessel independently. With this set-up both vessels can operate at different temperatures. This is the recommended kit.

- R/H upper return manifold with vent valve
- L/H upper return manifold with vent valve
- R/H lower supply manifold with drain valve
- · L/H lower supply manifold with drain valve
- 4 x quick-release insulated hoses, 60 cm
- · 4 x sealing caps for hoses



Single circulator manifold and hose kit (RR121305)

This single circulator kit is designed to work with one thermoregulator controlling the temperature of both vessels simultaneously. With this set-up both vessels would operate at the same temperature. The manifolds are supplied with blanking plugs should you wish to run only one vessel.

- · central upper return manifold with vent valve
- · central lower supply manifold with drain valve
- 4 x quick-release insulated hoses, 70 cm
- 4 x sealing caps for hoses
- 2 x stainless steel blanking plugs



Hei-Torque overhead stirrers

The powerful Hei-Torque range can accomplish the most demanding mixing tasks whilst providing the highest safety

Includes 3 YEAR warranty

375 375 2000

Stirrer models include

- Powerful stirring from 10 to 2000 rpm (Max. recommended speed with Reactor-Ready is 500 rpm).
- Available with basic (Value) or advanced (Precision) features, and a choice of up to 40, 200 or 400 Ncm torque.
- Intuitive control panel for easy operation, safe start and stop to avoid accidental start-up.
- New features for Precision range: set speed limit, adjustable acceleration, set a torque limit.
- New unique Quick-chuck for immediate and convenient 'one-hand' impellor changes, without tools.

	Hei-Torque Core	Hei-Torque Value 200	Hei-Torque Value 400	Hei-Torque Precision 200	Hei-Torque Precision 400
Power rating, motor input/output (W)	105/75	120/80	150/90	120/80	150/90
Number of gears	1	1	2	1	2
Speed range (rpm)	20 to 2000	10 to 2000	10 to 400 and 200 to 2000	10 to 2000	10 to 400 and 200 to 2000
Speed indicator	digital monochrome 2"	digital monochrome 2.4"	digital monochrome 2.4"	digital colour 3.2"	digital colour 3.2"
Speed control	electronic	electronic	electronic	electronic	electronic
Torque, maximum (Ncm)	40	200	400	200	400
Torque indicator (Ncm)	symbol	symbol	symbol	precise value	precise value
Overheat protection	automatic cut-out	automatic cut-out	automatic cut-out	automatic cut-out	automatic cut-out
Viscosity up to (mPa s)	10,000	100,000	250,000	100,000	250,000
Stirring cap. H ₂ O max. (L)	25	50	100	50	100
Shaft diameter, max. (mm)	10.5	10.5	10.5	10.5	10.5
Dimensions (W x H x D) (mm)	70 x 281.5 x 195	86 x 350 x 247	93 x 350 x 247	86 x 350 x 247	93 x 350 x 247
Weight (kg)	2.3	5.1	5.3	5.1	5.3
Set speed limit, Torque limit, Adjustable acceleration	-	-	-	Yes	Yes
Analogue/digital interface	-	-	-	USB and RS232	USB and RS232

For Reactor-Ready and Reactor-Ready Duo the maximum speed is 500 rpm for continuous use for 800 rpm for short periods.

Thermoregulators/circulators

Which thermoregulator should you choose?

- Radleys are a Huber distribution and service partner, so we recommend the Huber range of thermoregulators, circulators, chillers and waterbaths. However, Reactor-Ready is compatible with all leading brands of thermoregulator including Julabo, Lauda, Polyscience, Haake and many others.
- Also available are a range of hoses, adapters, temperature probes and thermofluid for all applications.

Description	Temperature range °C	Heating power (kW)	Cooling power at 0 °C (kW)	Application Examples		
Huber	45 °C to 200 °C	2.0	0.38	250 ml	Cool to -10 °C in <1 hr, then heat to 130 °C in <30 mins	
Ministat 230	at 230 -45 °C to 200 °C		0.30	1 Litre	Control at 0 °C, then heat to 100 °C in 20 mins	
Huber	-40 °C to 200 °C	1.5	0.45	500 ml	Cool to -20 °C in approx 30 mins, then heat to 150 °C in <30 mins	
Petite Fleur	-40 C to 200 C	1.0		1 Litre	Heat to 150 °C, then cool to 0 °C at > 2.5 K/min	
Huber	uber		0.7	1 Litre	Cool to -30 °C in <1 hr, then heat to 180 °C in <1 hr	
Unistat Tango -45 °C to 250 °C		3.0		2 Litre	Heat to 180 °C, then cool to 20 °C in approx 1 hr	
Huber	Huber 45 0 050 0		1.0	2 Litre	Cool to -20 °C in <30 mins, then heat to 100 °C in <30 mins	
Unistat 405	-45 °C to 250 °C	3.0	1.0	5 Litre	Control at 20 °C, then cool to -20 °C in approx 70 mins	
Huber 75 to 4 252 to		3.0	0.05	1 Litre	Control at 20 °C, then cool to -50 °C in 40 mins	
Unistat 705	-75 °C to 250 °C	3.0	0.65	5 Litre	Cool to -40 °C at approx 1 K/min, then heat to 150 °C in <1 hr	

The above applications are intended as a guide to the relative performance of each Huber thermoregulator with a range of vessel sizes and are not a prediction of actual performance. To achieve optimum performance the use of vacuum jacketed vessels is recommended for below ambient temperature applications. For comprehensive advice on thermoregulator selection please contact your local Radleys distributor.



Huber Unistat 405



Huber Petite Fleur - Baby Tango

Reaction vessels

Single jacketed, vacuum jacketed and process styles

Features include

- · Precision engineered borosilicate glass vessels can be used with Reactor-Ready and Reactor-Ready Duo.
- · DN100 ground and polished flat flange combines with the FEP-encapsulated O-ring and our unique self centering PTFE collar to offer a leak-tight seal between the vessel and lid.
- · Convenient and cost-effective vessel kits include stirrer, temperature probe and adapter.

Vessels feature zero dead space PTFE bottom outlet valve

- Exclusive stopcock design prevents over-tightening that can cause glass breakage.
- · Piston design allows quick and safe removal, eliminating the cleaning and assembly problems associated with other stopcock styles.
- Glass-filled PTFE piston minimises leaking or weeping across the full temperature range (-70 °C to +230 °C)
- Interchangeable between vessels with a full range of replacement parts and FEP O-rings.
- The large 15 mm bore minimises sample hold-up and is ideal for slurries (25 mm available on request).

Glass vessel lids

- DN100 5 neck glass lid fits all vessel sizes.
- DN100 lid with 5 x Rodaviss, grease-free, safety joints. B24 centre, B34, B29, B24 and B19 side joints.
- Rodaviss glass sockets will accept any glassware with traditional A, B or C length conical cones.

PTFE vessel lids

- DN100 7 neck PTFE lids fit all vessel sizes.
- · More robust and less likely to break than glass.
- Ideal for use with analytical probes.
- · Lids accept a range of removable Rodaviss sockets and compression fittings.
- · Rodaviss PTFE sockets will accept any glassware with traditional A, B or C length conical cones.

Standard Vessels

250 ml single iacketed vessel





1000 ml single iacketed vessel



250 ml vacuum Jacketed vessel



1000 ml vacuum iacketed vessel



500 ml single iacketed vessel



2000 ml single iacketed vessel



500 ml vacuum iacketed vessel



2000 ml vacuum iacketed vessel



5000 ml single iacketed vessel



5000 ml vacuum iacketed vessel



Process Vessels - mimic the geometry of plant scale reactors

Feature a 1.25 to 1 ratio of internal height to internal diameter, with dish-shaped bottom.

100 ml single jacketed process vessel



1000 ml single jacketed

process vessel

100 ml vacuum jacketed process vessel



1000 ml vacuum jacketed process vessel



250 ml single jacketed process vessel



2000 ml single jacketed process vessel



250 ml vacuum jacketed process vessel



2000 ml vacuum jacketed process vessel



500 ml single jacketed process vessel



5000 ml single jacketed process vessel



500 ml vacuum jacketed process vessel



5000 ml vacuum jacketed process vessel



Note: 3 litre vessels not shown

Reaction vessel kits

Reaction vessel kits make ordering vessels and accessories simple and much more cost effective







probe port B14

How to order Reactor-Ready

- 1. Select the Reactor-Ready core: RR121000.
- 2. Choose the vessel kit or kits you require by selecting the vessel volume and style.
- Select the overhead stirrer you require.
- If you need a thermoregulator, hoses, hose adapters, thermofluid or accessory glassware, then select from the accessory list opposite.
- 5. If you need automation, add the AVA Software and Data Hub.

Reactor-Ready core system

RR121000

Reactor-Ready core system includes:

- Stand with base, drip tray, support rods, bosses and all clamps
- Vessel to framework insulated hoses with quick-release connections
- Upper and lower drain-down manifolds
- · PTFE stirrer guide and drop-in coupling
- · Glass reaction vessel lid with Rodaviss sealing caps
- Vessel FEP O-ring seal and PTFE centering collar
- · 2 x Stainless steel vessel support collars

Important: Hoses to connect the Reactor-Ready framework to your thermoregulator are NOT included. You will require hoses with an M24x1.5 female connection or appropriate adapters. See opposite for hoses, adapters and suitable thermoregulators.

Standard Reaction Vessel Kits

Kits include: standard vessel, PTFE anchor stirrer paddle, PTFE Pt100 temp probe (with LEMO plug for connection to the thermoregulator) and PTFE temp probe adapter.

RR120010	250ml Single Jacketed Vessel Kit
RR120015	250ml Vacuum Jacketed Vessel Kit
RR120020	500ml Single Jacketed Vessel Kit
RR120025	500ml Vacuum Jacketed Vessel Kit
RR120030	1000ml Single Jacketed Vessel Kit
RR120035	1000ml Vacuum Jacketed Vessel Kit
RR120040	2000ml Single Jacketed Vessel Kit
RR120045	2000ml Vacuum Jacketed Vessel Kit
RR120042	3000ml Single Jacketed Vessel Kit
RR120047	3000ml Vacuum Jacketed Vessel Kit
RR120050	5000ml Single Jacketed Vessel Kit
RR120055	5000ml Vacuum Jacketed Vessel Kit



Kit includes Anchor Stirrer

Process Reaction Vessel Kits

Kits include: process vessel, PTFE turbine stirrer paddle, PTFE Pt100 temp probe (with LEMO plug for connection to the thermoregulator) and PTFE temp probe adapter.

RR120105 RR120107 RR120110	100ml Single Jacketed Process Vessel Kit 100ml Vacuum Jacketed Process Vessel Kit 250ml Single Jacketed Process Vessel Kit
RR120115	250ml Vacuum Jacketed Process Vessel Kit
RR120120	500ml Single Jacketed Process Vessel Kit
RR120125	500ml Vacuum Jacketed Process Vessel Kit
RR120130	1000ml Single Jacketed Process Vessel Kit
RR120135	1000ml Vacuum Jacketed Process Vessel Kit
RR120140	2000ml Single Jacketed Process Vessel Kit
RR120145	2000ml Vacuum Jacketed Process Vessel Kit
RR120142	3000ml Single Jacketed Process Vessel Kit
RR120147	3000ml Vacuum Jacketed Process Vessel Kit
RR120150	5000ml Single Jacketed Process Vessel Kit
RR120155	5000ml Vacuum Jacketed Process Vessel Kit



Kit includes **Turbine Stirrer**

How to order Reactor-Ready Duo

- 1. Select the Reactor-Ready Duo core: RR121300.
- 2. Choose manifold and hose kit for one or two circulators.
- Choose the vessel kit or kits you require by selecting the vessel volume and style.
- Select the overhead stirrers you require.
- If you need a thermoregulator, hoses, hose adapters, thermofluid or accessory glassware, then select from the accessory list opposite.
- 6. If you need automation, add the AVA Control Software and Data Hub.

Reactor-Ready Duo core system

Reactor-Ready Duo core system includes:

- Duo stand with base, drip tray, support rods, bosses and all clamps
- 2 x PTFE stirrer guides and drop-in couplings
- 2 x Glass reaction vessel lids with Rodaviss sealing caps
- 2 x Vessel FEP O-ring seal and PTFE centering collars
- 2 x Stainless steel vessel support collars

Important: Hoses to connect the Reactor-Ready Duo framework to your thermoregulator are NOT included. You will require hoses with an M24x1.5 female connection or appropriate adapters. See opposite for hoses, adapters and suitable thermoregulators.

Double Circulator Manifold and Hose Kit

Double Circulator Manifold and Hose Kit includes:

- R/H Upper Return Manifold with Vent Valve
- L/H Upper Return Manifold with Vent Valve
- R/H Lower Supply Manifold with Drain Valve
- L/H Lower Supply Manifold with Drain Valve
- 4 x Quick-Release Insulated Hoses, 60cm
- 4 x Sealing Caps for Hoses

Single Circulator Manifold and Hose Kit

Single Circulator Manifold and Hose Kit includes:

- Central Upper Return Manifold with Vent Valve
- Central Lower Supply Manifold with Drain Valve
- 4 x Quick-Release Insulated Hoses, 70cm
- 4 x Sealing Caps for Hoses
- · 4 x Stainless Steel Blanking Plugs

Hei-Torque Overhead Stirrers

501-60410-00 Hei-Torque Core 501-62010-00

Hei-Torque 200 Value 501-62020-00 Hei-Torque 200 Precision (USB/RS232)

501-64010-00 Hei Torque 400 Value

501-64020-00 Hei-Torque 400 Precision (USB/RS232)

General Maintenance Accessories

RR121150	Reactor-Ready Maintenance Kit (includes items listed below)
RR121027	Oil Wipe for Support Rods, Pack 10
RR121063	Stirrer Guide Bearing Seal B24 10mm
RR121064	Stirrer Guide Ferrule and Shaft Seal B24 10mm
RR121071	23mm Viton O-Rings for Quick Release Hose Coupling, Pack 2
RR121072	28mm Viton O-Rings for Quick Release Hose Coupling, Pack 2
RR121090	FEP Encapsulated Silicone Vessel O-Ring
RR166100	Replacement 15mm Piston Top O-Ring
RR166102	Replacement 15mm Piston Middle O-Ring

Reactor-Ready Duo reaction vessel compatibility guide

The following combinations of reaction vessel sizes can be used in the Reactor-Ready Duo.

	100ml	250 ml	500 ml	1000 ml	2000 ml	3000 ml	5000 ml	
Standard Vessel, Single Jacketed								
250 ml	n/a	•	•	•	•	•	•	
500 ml	n/a	•	•	•	•	•	•	
1000 ml	n/a	•	•	•	•	•	•	
2000 ml	n/a	•	•	•	•	•	•	
5000 ml	n/a	•	•	•	•	•	•	
		Proce	ss Vessel,	Single Jac	keted			
100 ml	•	•	•	•	•	•	•	
250 ml	•	•	•	•	•	•	•	
500 ml	•	•	•	•	•	•	•	
1000 ml	•	•	•	•	•	•	•	
2000 ml	•	•	•	•	•	•	•	
3000 ml	•	•	•	•	•	•	•	
5000 ml	•	•	•	•	•	n/a	n/a	

	100ml	250 ml	500 ml	1000 ml	2000 ml	3000 ml	5000 ml	
Standard Vessel, Vacuum Jacketed								
250 ml	n/a	•	•	•	•	•	•	
500 ml	n/a	•	•	•	•	•	•	
1000 ml	n/a	•	•	•	•	•	•	
2000 ml	n/a	•	•	•	•	•	•	
5000 ml	n/a	•	•	•	•	n/a	n/a	
		Proces	ss Vessel, \	Vacuum Ja	cketed			
100 ml	•	•	•	•	•	•	•	
250 ml	•	•	•	•	•	•	•	
500 ml	•	•	•	•	•	•	•	
1000 ml	•	•	•	•	•	n/a	n/a	
2000 ml	•	•	•	•	•	n/a	n/a	
3000 ml	•	•	•	•	•	•	•	
5000 ml	•	•	•	n/a	n/a	n/a	n/a	

```
Glass DN100 Vessel Lid
Huber Thermoregulators (Circulators)
                                                                                                                                                        5 Neck Glass Lid - 1 x B19, 2 x B24, 1 x B29, 1 x B34
                             Ministat 240 Pilot ONE Controller
HB2016 0005 01
                                                                                                                                   RR136000
                                                                                                                                                        Includes Rodaviss sealing caps
HR1030 0001 01
                              Petite Fleur Pilot ONE Controller
HB1030 0001 01FT
                             Petite Fleur with 4 Litre Expansion Tank
                                                                                                                                   PTFE DN100 Vessel Lids
                             Unistat Tango Pilot ONE Controller
HB1000 0037 01
                                                                                                                                                        5 Neck PTFE Lid with 2 Blanked Ports. 1 x B19, 2 x B24, 1 x B29, 1 x B34 Includes Rodaviss sealing caps
                                                                                                                                   RR136100
                             Unistat 405 Pilot ONE Controller
HB1002.0045.01
                             Unistat 705 Pilot ONE Controller (230V or 400V)
HB1001.0041.01
                                                                                                                                                        7 Neck PTFE Lid with 7 Threaded Ports. 2 x 1/4" NPT, 5 x 3/4" NPT PTFE inserts not included - select from below
                                                                                                                                   RR136110
Thermoregulator Hoses, Adapters and Thermofluid
HR9325
                     M24x1.5 Insulated Hose NW12 100cm -60°C to +260°C
                                                                                                                                   PTFE Rodaviss Sockets, Compression Fittings and Blanking Plugs for PTFE Lids
                     M24x1.5 Insulated Hose NW12 150cm -60°C to +260°C M24x1.5 Insulated Hose NW12 200cm -60°C to +260°C -60°C to +260°C
                                                                                                                                                        Rodaviss B19 Socket Adapter (
Rodaviss B24 Socket Adapter (
HR9326
HB9327
                                                                                                                                                                                                     NIPT
                                                                                                                                   RR136301
                                                                                                                                   RR136302
                                                                                                                                                         Rodaviss B29 Socket Adapter
HB9328
                     M24x1.5 Insulated Hose NW12 300cm -60°C to +260°C
                                                                                                                                   RR136303
                                                                                                                                                        Rodaviss B34 Socket Adapter (3/4" NPT Rodaviss B45 Socket Adapter (3/4" NPT
Insulated Hose with Metal Inner
                                                                                                                                   RR136304
HR9274
                      M24x1.5 Insulated Hose NW12 100cm -100°C to +350°C
                                                                                                                                   RR136314
                                                                                                                                                         3mm Bore Compression Fitting 1/4" NPT
                      M24x1.5 Insulated Hose NW12 150cm -100°C to +350°C
HB9275
                                                                                                                                                        6mm Bore Compression Fitting ½" NPT 8mm Bore Compression Fitting ½" NPT
                                                                                                                                   RR136315
HB9276
                      M24x1.5 Insulated Hose NW12 200cm -100°C to +350°C
                                                                                                                                   RR136316
                                                                                                                                                        12mm Bore Compression Fitting 3/4 NPT 14mm Bore Compression Fitting 3/4 NPT 14mm Bore Compression Fitting 3/4 NPT
HB9277
                     M24x1.5 Insulated Hose NW12 300cm -100°C to +350°C
                                                                                                                                   RR136317
                                                                                                                                   RR136318
HB6724
                      M24x1.5 Female to M16x1 Male Adapter
                                                                                                                                   RR136319
                                                                                                                                                         16mm Bore Compression Fitting 3/4" NP
                                                                                                                                   RR136320
RR136327
                                                                                                                                                        19mm Bore Compression Fitting 3/4" NPT PTFE Blanking Plug 1/4" NPT
HR6945
                     M24x1.5 Male to M16x1 Female Adapter
HB6723
                      M24x1.5 Female to M30x1.5 Male Adapter
                                                                                                                                   RR136328
                                                                                                                                                        PTFE Blanking Plug 3/4" NPT
HB6195
                      M16x1 90° Adapter
HB9256
                      M24x1.5 90° Adapter
                                                                                                                                   Standard Reaction Vessels Only
HB6089
                      M16x1 Nut
                                                                                                                                   RR130250
RR130255
                                                                                                                                                        250ml Single Jacketed Vessel
250ml Vacuum Jacketed Vessel
                      M24x1.5 Nut
HB12634
HB6086
                      M16x1 Hose Connector NW8
                                                                                                                                   RR130500
                                                                                                                                                         500ml Single Jacketed Vessel
HB9471
                      M24x1.5 Hose Connector NW12
                                                                                                                                   RR130550
                                                                                                                                                        500ml Vacuum Jacketed Vessel
1000ml Single Jacketed Vessel
                     M16x1 Isolating Ball Valve -20°C...+140°C (max. 6 bar at +140°C) M24x1.5 Isolating Ball Valve -10°C...+180°C (max. 6 bar at +180°C)
HB6091
                                                                                                                                   RR131000
                                                                                                                                   RR131050
RR132000
                                                                                                                                                        1000ml Vacuum Jacketed Vessel
2000ml Single Jacketed Vessel
HB9236
 Thermofluid
                                                                                                                                    RR132050
                                                                                                                                                         2000ml Vacuum Jacketed Vessel
HB6161
                      SilOil M20.195/235.20. 5 Litre
                                                                   -20°C to +195°C/235°C
                                                                                                                                                        3000ml Single Jacketed Vessel
3000ml Vacuum Jacketed Vessel
                                                                 -20°C to +195°C/235°C
-40°C to +165°C/220°C
                                                                                                                                   RR133000
                     SilOil M20.195/235.20, 10 Litre
SilOil M40.165/220.10, 5 Litre
HR6162
                                                                                                                                   RR133050
HR6163
                                                                                                                                    RR135000
                                                                                                                                                         5000ml Single Jacketed Vessel
                                                                  -40°C to +165°C/220°C
HB6164
                      SilOil M40.165/220.10, 10 Litre
                                                                                                                                   RR135050
                                                                                                                                                        5000ml Vacuum Jacketed Vessel
                                                                  -60°C to +115°C/200°C
HB6165
                      SilOil M60.115/200.05. 5 Litre
                                                                                                                                  PTFE Stirrer Shafts (for Standard Vessels)
RR150250 Anchor PTFE Stirrer Shaft - 250ml vessels
RR150260 Retreat Curve PTFE Stirrer Shaft - 250ml vessels
Turbine PTFE Stirrer Shaft - 250ml vessels
                                                                 -60°C to +115°C/200°C
HB6166
                      SilOil M60.115/200.05, 10 Litre
HB6479
                      DW-Therm M90.200.02, 10 Litre -90°C to +200°C (Unistats only)
Glassware Accessories
                                                                                                                                                        Anchor PTFE Stirrer Shaft - 500ml vessels
Retreat Curve PTFE Stirrer Shaft - 500ml vessels
Turbine PTFE Stirrer Shaft - 500ml vessels
RR139005
                     Condenser Liebig 250mm B29 + GL14 + fittings
                                                                                                                                   RR150500
RR139007
                     Condenser Double Surface 200mm B29 + GL14 + fittings
                                                                                                                                   RR150505
                                                                                                                                   RR150510
RR139009
                     Condenser Jacket Coil Rodaviss B29 + GL14 + fittings
                                                                                                                                                        Anchor PTFE Stirrer Shaft - 1000ml vessels
Retreat Curve PTFE Stirrer Shaft - 1000ml vessels
Turbine PTFE Stirrer Shaft - 1000ml vessels
                                                                                                                                   RR151000
RR139011
                     Right Angle Adapter B29 + GL14 + PTFE fittings
                                                                                                                                   RR151005
                     50ml Graduated Dropping Funnel Equalising Arm B24
RR139013
                                                                                                                                   RR151010
                     100ml Graduated Dropping Funnel Equalising Arm B24
250ml Graduated Dropping Funnel Equalising Arm B24
RR139015
                                                                                                                                                        Anchor PTFE Stirrer Shaft - 2000ml vessels
Retreat Curve PTFE Stirrer Shaft - 2000ml vessels
Turbine PTFE Stirrer Shaft - 2000ml vessels
                                                                                                                                   RR152000
RR139017
                                                                                                                                   RR152005
                     500ml Graduated Dropping Funnel Equalising Arm B24
1000ml Graduated Dropping Funnel Equalising Arm B24
RR139019
                                                                                                                                   RR152010
                                                                                                                                                        Anchor PTFE Stirrer Shaft - 3000ml vessels
Retreat Curve PTFE Stirrer Shaft - 3000ml vessels
Turbine PTFE Stirrer Shaft - 3000ml vessels
Anchor PTFE Stirrer Shaft - 3000ml vessels
Anchor PTFE Stirrer Shaft - 5000ml vessels
Retreat Curve PTFE Stirrer Shaft - 5000ml vessels
Turbine PTFE Stirrer Shaft - 5000ml vessels
RR139021
                                                                                                                                   RR153000
RR153005
                     Swan Neck Adapter Angled B24 to Vertical B24
RR139023
RR139025
                     Gas Purge Adapter B24 + GL14 + fittings
                                                                                                                                   RR153010
                                                                                                                                   RR155000
RR139027
                     Glass Cone Stopper B19
RR139029
                     Glass Cone Stopper B24
                                                                                                                                   RR155010
                                                                                                                                                         Turbine PTFE Stirrer Shaft - 5000ml vessels
RR139031
                     Glass Cone Stopper B29
RR139033
                     Glass Cone Stopper B34
                                                                                                                                   Process Reaction Vessels Only
                     Enlarging Adapter B24 to B29
Enlarging Adapter B29 to B34
Reducing Adapter B24 to B19
                                                                                                                                                         100ml Single Jacketed Process Vessel
100ml Vacuum Jacketed Process Vessel
RR139035
                                                                                                                                   RR130102
RR130152
RR139037
RR139039
                                                                                                                                   RR130252
RR130257
                                                                                                                                                        250ml Single Jacketed Process Vessel 250ml Vacuum Jacketed Process Vessel
                     Reducing Adapter B29 to B19
Reducing Adapter B29 to B24
RR139041
                                                                                                                                                        500ml Single Jacketed Process Vessel
500ml Vacuum Jacketed Process Vessel
1000ml Single Jacketed Process Vessel
RR139043
                                                                                                                                   RR130502
RR130552
RR139045
                     Reducing Adapter B34 to B19
RR139047
                     Reducing Adapter B34 to B24
                                                                                                                                   RR131052
RR132002
                                                                                                                                                        1000ml Vacuum Jacketed Process Vessel
2000ml Single Jacketed Process Vessel
RR139049
                     Reducing Adapter B34 to B29
RR139051
                     Oil Bubbler Type 1 B24 + GL14 + fittings
                                                                                                                                   RR132052
                                                                                                                                                         2000ml Vacuum Jacketed Process Vessel
RR139053
                     Solid Additions Funnel B34 (with extended tube)
                                                                                                                                                        3000ml Single Jacketed Process Vessel
3000ml Vacuum Jacketed Process Vessel
                                                                                                                                   RR133002
RR139057
                     Valved Liquid Inlet Adapter B24 to GL14
                                                                                                                                   RR133052
RR139063
                     PTFE Hose Adapter B19 for 12mm ID Tubing
                                                                                                                                   RR135002
                                                                                                                                                        5000ml Single Jacketed Process Vessel
                                                                                                                                   RR135052
RR139067
                                                                                                                                                        5000ml Vacuum Jacketed Process Vessel
                     Boss Head
                     3-Prong Clamp
PTFE Flexible Bellows B24
RR139069
                                                                                                                                   PTFE Stirrer Shafts (for Process Vessels)
RR150102 Anchor PTFE Stirrer Shaft - 100ml Process Vessels
RR150112 Retreat Curve PTFE Stirrer Shaft - 100ml Process Vessels
Turbine PTFE Stirrer Shaft - 100ml Process Vessels
Turbine PTFE Stirrer Shaft - 100ml Process Vessels

RR150112 STEE STIRRER STAFT - 100ml Process Vessels

RR150112 STEE STIRRER STAFT - 100ml Process Vessels
RR139079
                     PTFE Flexible Bellows B29
RR139081
                     PTFE Flexible Bellows B34
RR139083
RR139150
                     Reflux Divider Head B29 + GL18 + PTFE Fitting
                                                                                                                                                        Anchor PTFE Stirrer Shaft - 250ml Process Vessels
Retreat Curve PTFE Stirrer Shaft - 250ml Process Vessels
Turbine PTFE Stirrer Shaft - 250ml Process Vessels
                                                                                                                                    RR150252
RR139152
                     Reflux Divider Head B29 + GL18 + Temp Probe Port B14
                                                                                                                                   RR150257
RR150262
RR139190
                     Reflux Divider Kit
                                                                                                                                                        Anchor PTFE Stirrer Shaft - 500ml Process Vessels
Retreat Curve PTFE Stirrer Shaft - 500ml Process Vessels
Turbine PTFE Stirrer Shaft - 500ml Process Vessels
RR139192
                     Reflux Divider Kit with Vapour Temp Kit
                                                                                                                                   RR150502
RR150507
                     250ml Receiver Flask 2 \dot{x} B24 Socket, + B24 Cone and Tap
RR139154
                                                                                                                                   RR150512
                     500ml Receiver Flask 3 x B24 Socket, + B24 Cone and Tap
RR139324
                                                                                                                                                        Anchor PTFE Stirrer Shaft - 1000ml Process Vessels
Retreat Curve PTFE Stirrer Shaft - 1000ml Process Vessels
                                                                                                                                   RR151002
                     1000ml Receiver Flask 3 x B24 Socket, + B24 Cone and Tap
RR139326
                                                                                                                                   RR151007
                                                                                                                                                        Turbine PTFE Stirrer Shaft - 1000ml Process Vessels
Anchor PTFE Stirrer Shaft - 2000ml Process Vessels
Retreat Curve PTFE Stirrer Shaft - 2000ml Process Vessels
                                                                                                                                   RR151012
PTFE Pt100 Temperature Probes and Adapters
                     Probe 8mm with LEMO Plug - 100ml, 250ml + 500ml Vessels
Probe 9mm with LEMO Plug - 1000ml, 2000ml + 3000ml Vessels
                                                                                                                                   RR152002
RR140250
                                                                                                                                    RR152007
RR141000
                                                                                                                                                         Turbine PTFE Stirrer Shaft - 2000ml Process Vessels
Anchor PTFE Stirrer Shaft - 3000ml Process Vessels
                                                                                                                                   RR152012
RR145000
                     Probe 9mm with LEMO Plug - 5000ml Vessel
                                                                                                                                   RR153002
RR140010
                     PTFE Temperature Probe Adapter - 8mm
                                                                                                                                   RR153007
                                                                                                                                                         Retreat Curve PTFE Stirrer Shaft - 3000ml Process Vessels
RR140015
                     PTFE Temperature Probe Adapter - 9mm
                                                                                                                                                        Turbine PTFE Stirrer Shaft - 3000ml Process Vessels
Anchor PTFE Stirrer Shaft - 5000ml Process Vessels
Anchor PTFE Stirrer Shaft - 5000ml Process Vessels
Retreat Curve PTFE Stirrer Shaft - 5000ml Process Vessels
Turbine PTFE Stirrer Shaft - 5000ml Process Vessels
                                                                                                                                   RR155002
Reaction Vessel Support Collar
                                                                                                                                   RR155007
                                                                                                                                   RR155012
RR121048
                     Vessel Support Collar
```

AVA Lab Control Software - take control of your chemistry

Control multi-device jacketed reaction system

Multi-device reactions

- Control up to 4 reaction systems with up to 16 devices on one screen
- · Create complex experiments with any number of steps in series or in parallel
- · Pre-program multi-step recipe, with the flexibility to make and track on-the-fly adjustments
- · Interlink devices and set feedback/control loops, end-point conditions and safety limits

No more manual data logging or manual errors

- · Automatically record reaction parameters and log what you do, as you do it
- · Repeat experiments accurately for reproducible and consistent results
- · Create reports in a few clicks or export data as a CSV file for further analysis
- · Share results between users to improve research and collaboration

Safe unattended chemistry

- · Automate cooling during exothermic events
- · Link devices such as balances and pumps for controlled, unattended reagent additions
- · Define safety overrides and cut-off conditions
- · Configure audio and visual alarms

Data Hub

Using the Data Hub, AVA can control and log multiple 3rd party devices with RS232 serial interface.



Data Hub Specifications

- 4 x RS232 ports
- 2 x Pt100 temperature sensor ports with LEMO connection
- 1 x Ethernet port
- 2 x Ethernet cables

AVA Software Kit includes:

- AVA Software
- Laptop controller
- 4 Port data hub
- 1 year AVA Care Support

AVA Care Support

- · Free support for 1st year
- · Free priority email and telephone support
- · Free set up and application support
- Free software updates during support period, keeping software current

Download and try the software for FREE

Learn how AVA software works - try before you buy

- · Simulate control of devices
- · Set up apparatus and control experimental Schedules
- · Share setups and schedules with other AVA users
- · Analyse results and create reports for real or simulated data
- Find out more about AVA Level 1-3 on our website

Download Free AVA Demo www.radleys.com







Control a variety of devices

AVA software includes a library of pre-configured driver files allowing easy integration with a wide range of 3rd party devices.

Radleys also provide a 'New Driver Configuration Service' if required.
























Radleys provides innovative chemistry equipment for safer, cleaner, greener and more productive chemical research.

We have been manufacturing scientific glassware and laboratory instruments for over 50 years and our customers include leading blue-chip industrial and academic research facilities around the world. Our areas of expertise are focused on equipment for chemical synthesis, process development, work-up and evaporation.

Visit www.radleys.com to see our full range of chemistry productivity tools.



Shire Hill, Saffron Walden, Essex, CB11 3AZ. United Kingdom. t: +44 1799 513320 f: +44 1799 513283 e: sales@radleys.co.uk