

### BAKERBOND™ Chiral Columns

BAKERBOND™ Chiral Phase HPLC columns offer economic alternatives for the separation of racemic mixtures, including isomers of pharmaceuticals and other biologically active compounds.

#### CHIRAL HPLC COLUMNS—PIRKLE-TYPE

Description	Size	Column Size	Particle Shape	Product Number
® DNBP <sup>1</sup> (Ionic)	4.6 x 250 mm	5 µm	Spherical	7103-00
® DNBP <sup>1</sup> (Covalent)	4.6 x 250 mm	5 µm	Spherical	7113-00
	21.2 x 250 mm	10 µm	Irregular	7153-00
(S) DNBL <sup>2</sup> (Covalent)	4.6 x 250 mm	5 µm	Spherical	7115-00
α-Naphthyl Urea	4.6 x 250 mm	5 µm	Spherical	7309-00

<sup>1</sup> ®-N-3,5-Dinitrobenzoyl-phenyl glycine

<sup>2</sup> (S)-N-3,5-Dinitrobenzoyl-leucine

#### CHIRAL HPLC COLUMNS—PROTEIN-TYPE

Description	Size	Column Size	Particle Shape	Product Number
CHIRAL-AGP™	3 x 10 mm	5 µm	Spherical	7400-00
	4 x 50 mm	5 µm	Spherical	7165-25
	4 x 100 mm	5 µm	Spherical	7165-00
	3 x 100 mm	5 µm	Spherical	7165-01
	3 x 150 mm	5 µm	Spherical	7165-02

**tyco** / Healthcare / **Mallinckrodt Baker**

www.jtbaker.com

<b>Mallinckrodt Baker</b> <b>Holland Headquarters</b> Teugseweg 20 PO Box 1 7400 AA Deventer The Netherlands t 0031 570 687500 f 0031 570 687574 e jtbaker.nl@ emea.tycohealthcare.com	<b>Mallinckrodt Baker</b> <b>UK</b> 107/112 Leaderhall Street London EC3A 4AH United Kingdom t 0044 1908 506000 f 0044 1908 503290 e jtbaker.uk@ emea.tycohealthcare.com	<b>Mallinckrodt Baker</b> <b>Italy</b> Via P. Portoluppi 11/2 20138 Milano Italy t 0039 02 580891 f 0039 02 58019264 e jtbaker.it@ emea.tycohealthcare.com	<b>Mallinckrodt Baker</b> <b>Germany</b> Im Leuschnerpark 4 D-64347 Griesheim Germany t 0049 6155 667939 f 0049 6155 667940 e jtbaker.de@ emea.tycohealthcare.com	<b>Mallinckrodt Baker</b> <b>France</b> 33 Avenue Philippe Auguste 75011 Paris France t 0033 144 64 89 71/2 f 0033 144 64 89 73 e jtbaker.fr@ emea.tycohealthcare.com
---	---	--	---	---



**A CENTURY  
OF EMPOWERING  
DISCOVERY™**

### CHIRAL HPLC COLUMNS—PROTEIN-TYPE – continued

Description	Size	Column Size	Particle Shape	Product Number
CHIRAL-AGP™	2 x 100 mm	5 µm	Spherical	7165-03
	2 x 150 mm	5 µm	Spherical	7165-04
	4 x 50 mm	5 µm	Spherical	7165-26
	4 x 150 mm	5 µm	Spherical	7165-30
	10 x 100 mm	5 µm	Spherical	7404-37
	10 x 150 mm	5 µm	Spherical	7404-38
CHIRAL-HAS™	3 x 10 mm	5 µm	Spherical	7401-00
	3 x 100 mm	5 µm	Spherical	7401-01
	4 x 100 mm	5 µm	Spherical	7401-29
	4 x 150 mm	5 µm	Spherical	7401-30
Guard Column Holder	—	—	—	7402-00